

A5158

Leader in Biomolecular Solutions for Life Science



SUN5 Rabbit pAb

Catalog No.: A5158

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa/42kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to enable protein-membrane adaptor activity. Involved in spermatid development. Located in sperm connecting piece. Human ortholog(s) of this gene implicated in spermatogenic failure 16. Orthologous to human SUN5 (Sad1 and UNC84 domain containing 5).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

76407

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SRG4; Spag4l; SPAG4L-2; 1700021O15Rik; SUN5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.