

A5178

Leader in Biomolecular Solutions for Life Science



# JAKMIP1 Rabbit pAb

Catalog No.: A5178

## Basic Information

---

**Observed MW**

85kDa

**Calculated MW**

73kDa

**Category**

Polyclonal Antibody

**Applications**

WB

**Cross-Reactivity**

Human,Mouse,Rat

## Background

---

Enables GABA receptor binding activity and RNA binding activity. Involved in cognition. Is extrinsic component of membrane. Part of ribonucleoprotein complex.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

**Gene ID**

152789

**Swiss Prot**

Q96N16

**Immunogen**

Synthetic peptide. This information is considered to be commercially sensitive.

**Synonyms**

JAMIP1; MARLIN1; Gababrbp; JAKMIP1

## Contact

---

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

**Source**

Rabbit

**Isotype**

IgG


**Purification**

Affinity purification

**Storage**

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,pH7.3.

## Validation Data



Western blot analysis of various lysates using JAKMIP1 Rabbit pAb (A5178) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.