

A5208

Leader in Biomolecular Solutions for Life Science



TMED4 Rabbit pAb

Catalog No.: A5208

Basic Information

Observed MW

Refer to figures

Calculated MW

26kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Involved in positive regulation of I-kappaB kinase/NF-kappaB signaling. Predicted to be located in endoplasmic reticulum membrane. Predicted to be integral component of membrane. Predicted to be active in several cellular components, including COPII-coated ER to Golgi transport vesicle; Golgi apparatus; and endoplasmic reticulum-Golgi intermediate compartment.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

222068

Swiss Prot

Q7Z7H5

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HNLF; ERS25; p24a3; GMP25iso; p24alpha3; TMED4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.