

A5575

Leader in Biomolecular Solutions for Life Science



RABL2A Guinea pig pAb

Catalog No.: A5575

Basic Information

Observed MW

Refer to figures

Calculated MW

26kDa/27kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables GTP binding activity. Involved in photoreceptor cell morphogenesis. Acts upstream of or within cilium assembly; flagellated sperm motility; and single fertilization. Located in sperm midpiece. Part of intraciliary transport particle B. Is expressed in several structures, including eye; genitourinary system; gut; nervous system; and respiratory system. Orthologous to several human genes including RABL2B (RAB, member of RAS oncogene family like 2B).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

68708

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Rabl2a; 1110031N17Rik

Contact



www.abclonal.com

Product Information

Source

Guinea pig

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.