

A5721

Leader in Biomolecular Solutions for Life Science



FUT2 Rabbit pAb

Catalog No.: A5721

Basic Information

Observed MW

35kDa

Calculated MW

39kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene is one of two encoding the galactoside 2-L-fucosyltransferase enzyme. The encoded protein is important for the final step in the soluble ABO blood group antigen synthesis pathway. It is also involved in cell-cell interaction, cell surface expression, and cell proliferation. Mutations in this gene are a cause of the H-Bombay blood group where red blood cells lack the H antigen.

Recommended Dilutions

WB 1:500 - 1:1000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

2524

Swiss Prot

Q10981

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SE; Se2; sej; SEC2; B12QTL1; FUT2

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of lysates from HT-29 cells, using FUT2 Rabbit pAb (A5721) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.