

A6275

Leader in Biomolecular Solutions for Life Science



# LRP12 Rabbit pAb

Catalog No.: A6275

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

95kDa

### Category

Polyclonal Antibody

### Applications

WB,IF/ICC,FC

### Cross-Reactivity

Human

## Background

This gene encodes a member of the low-density lipoprotein receptor related protein family. The product of this gene is a transmembrane protein that is differentially expressed in many cancer cells. Alternate splicing results in multiple transcript variants.

## Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:100
FC	1:20 - 1:50

## Immunogen Information

**Gene ID**  
29967

**Swiss Prot**  
Q9Y561

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

ST7; MIG13A; LRP12

## Contact

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

**Source**  
Rabbit

**Isotype**  
IgG

**Purification**  
Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide,pH7.3.