

A6448

Leader in Biomolecular Solutions for Life Science



## C1D Rabbit pAb

Catalog No.: A6448

### Basic Information

#### Observed MW

Refer to Figures

#### Calculated MW

16kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

The protein encoded by this gene is a DNA binding and apoptosis-inducing protein and is localized in the nucleus. It is also a Rac3-interacting protein which acts as a corepressor for the thyroid hormone receptor. This protein is thought to regulate TRAX/Translin complex formation. Alternate splicing results in multiple transcript variants that encode the same protein. Multiple pseudogenes of this gene are found on chromosome 10.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

10438

#### Swiss Prot

Q13901

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-141 of human C1D (NP\_001177194.1).

#### Synonyms

LRP1; hC1D; Rrp47; SUNCOR; SUN-CoR; C1D

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.