

A6484

Leader in Biomolecular Solutions for Life Science



POMGNT2 Rabbit pAb

Catalog No.: A6484

Basic Information

Observed MW

Calculated MW

67kDa

Category

Polyclonal Antibody

Applications

ELISA

Cross-Reactivity

Human

Background

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

Recommended Dilutions

Immunogen Information

Gene ID

84892

Swiss Prot

Q8NAT1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 301-580 of human POMGNT2 (NP_116195.2).

Synonyms

AGO61; GTDC2; MDDGA8; MDDGC8; C3orf39; POMGNT2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.