

A6494

Leader in Biomolecular Solutions for Life Science



## ZBTB38 Rabbit pAb

Catalog No.: A6494

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

134kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

The protein encoded by this gene is a zinc finger transcriptional activator that binds methylated DNA. The encoded protein can form homodimers or heterodimers through the zinc finger domains. In mouse, inhibition of this protein has been associated with apoptosis in some cell types.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

253461

#### Swiss Prot

Q8NAP3

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 896-1195 of human ZBTB38 (NP\_001073881.2).

#### Synonyms

CIBZ; ZNF921; PPP1R171; ZBTB38

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.