

A6592

Leader in Biomolecular Solutions for Life Science



# FRZB Rabbit pAb

Catalog No.: A6592

## Basic Information

---

### Observed MW

36kDa

### Calculated MW

36kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Rat,Mouse

## Background

---

The protein encoded by this gene is a secreted protein that is involved in the regulation of bone development. Defects in this gene are a cause of female-specific osteoarthritis (OA) susceptibility.

## Recommended Dilutions

---

WB 1:2000 - 1:7000

## Immunogen Information

---

### Gene ID

2487

### Swiss Prot

Q92765

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 166-325 of human FRZB (NP\_001454.2).

### Synonyms

FRE; OS1; FZRB; hFIZ; FRITZ; FRP-3; FRZB1; SFRP3; SRFP3; FRZB-1; FRZB-PEN; FRZB

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of lysates from Rat spleen, using FRZB Rabbit pAb (A6592) at 1:6000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25 $\mu$ g per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.