

A6860

Leader in Biomolecular Solutions for Life Science



# SREBP2 Rabbit pAb

Catalog No.: A6860

## Basic Information

### Observed MW

53kDa

### Calculated MW

124kDa

### Category

Polyclonal Antibody

### Applications

### Cross-Reactivity

Human,Mouse,Rat

## Background

This gene encodes a member of the a ubiquitously expressed transcription factor that controls cholesterol homeostasis by regulating transcription of sterol-regulated genes. The encoded protein contains a basic helix-loop-helix-leucine zipper (bHLH-Zip) domain and binds the sterol regulatory element 1 motif. Alternate splicing results in multiple transcript variants.

## Recommended Dilutions

## Immunogen Information

### Gene ID

6721

### Swiss Prot

Q12772

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

SREBP2; bHLHd2; SREBP-2

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.