

A6863

Leader in Biomolecular Solutions for Life Science



## AP2A1 Rabbit pAb

Catalog No.: A6863

1 Publications

### Basic Information

#### Observed MW

108kDa

#### Calculated MW

108kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse, Rat

### Background

This gene encodes the alpha 1 adaptin subunit of the adaptor protein 2 (AP-2) complex found in clathrin coated vesicles. The AP-2 complex is a heterotetramer consisting of two large adaptins (alpha or beta), a medium adaptin (mu), and a small adaptin (sigma). The complex is part of the protein coat on the cytoplasmic face of coated vesicles which links clathrin to receptors in vesicles. Alternative splicing of this gene results in two transcript variants encoding two different isoforms. A third transcript variant has been described, but its full length nature has not been determined.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

160

#### Swiss Prot

O95782

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human AP2A1 (NP\_055018.2).

#### Synonyms

ADTAA; CLAPA1; AP2-ALPHA; AP2A1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

# Validation Data



Western blot analysis of extracts of various cell lines, using AP2A1 antibody (A6863) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 10s.