

A7384

Leader in Biomolecular Solutions for Life Science



FASTKD1 Rabbit pAb

Catalog No.: A7384

Basic Information

Observed MW

110kDa

Calculated MW

97kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables RNA binding activity. Involved in mitochondrial RNA metabolic process and regulation of mitochondrial mRNA stability. Located in mitochondrion.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79675

Swiss Prot

Q53R41

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 638-847 of human FASTKD1 (NP_078898.3).

Synonyms

FASTKD1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using FASTKD1 antibody (A7384) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.