

A7583

Leader in Biomolecular Solutions for Life Science



WDR47 Rabbit pAb

Catalog No.: A7583

Basic Information

Observed MW

Refer to figures

Calculated MW

102kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to be located in cytoplasm and microtubule.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

22911

Swiss Prot

O94967

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

WDR47

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.