

A7944

Leader in Biomolecular Solutions for Life Science



SSNA1 Rabbit pAb

Catalog No.: A7944

Basic Information

Observed MW

14kDa

Calculated MW

14kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Enables identical protein binding activity. Predicted to act upstream of or within ciliary receptor clustering involved in smoothened signaling pathway and intraciliary transport. Located in centrosome.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

8636

Swiss Prot

O43805

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

N14; NA14; NA-14; SSNA1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.