

A7987

Leader in Biomolecular Solutions for Life Science



# CALCOCO1 Rabbit pAb

Catalog No.: A7987

1 Publications

## Basic Information

### Observed MW

77kDa

### Calculated MW

77kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse

## Background

Enables several functions, including armadillo repeat domain binding activity; beta-catenin binding activity; and nuclear receptor coactivator activity. Involved in positive regulation of gene expression and positive regulation of transcription, DNA-templated. Located in cytosol and nucleus.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

57658

### Swiss Prot

Q9P1Z2

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 462-691 of human CALCOCO1 (NP\_065949.1).

### Synonyms

Cocoa; PP13275; calphoglin; CALCOCO1

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data



Western blot analysis of extracts of various cell lines, using CALCOCO1 antibody (A7987) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.