

A8151

Leader in Biomolecular Solutions for Life Science



# SRM Rabbit pAb

Catalog No.: A8151

## Basic Information

### Observed MW

35kDa

### Calculated MW

34kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse,Rat

## Background

The polyamines putrescine, spermine, and spermidine are ubiquitous polycationic mediators of cell growth and differentiation. Spermidine synthase is one of four enzymes in the polyamine-biosynthetic pathway and carries out the final step of spermidine biosynthesis. This enzyme catalyzes the conversion of putrescine to spermidine using decarboxylated S-adenosylmethionine as the cofactor.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

6723

### Swiss Prot

P19623

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human SRM (NP\_003123.2).

### Synonyms

PAPT; SPS1; SPDSY; SRML1; SRM

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data



Western blot analysis of extracts of various cell lines, using SRM antibody (A8151) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 10s.