

A8234

Leader in Biomolecular Solutions for Life Science



KANSL3 Rabbit pAb

Catalog No.: A8234 **1 Publications**

Basic Information

Observed MW

120kDa

Calculated MW

96kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Involved in histone H4-K16 acetylation; histone H4-K5 acetylation; and histone H4-K8 acetylation. Located in nucleoplasm. Part of histone acetyltransferase complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

55683

Swiss Prot

Q9P2N6

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

NSL3; Rcd1; KIAA1310; KANSL3

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of various lysates using KANSL3 Rabbit pAb (A8234) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 15s.