

A8234

Leader in Biomolecular Solutions for Life Science



KANSL3 Rabbit pAb

Catalog No.: A8234 **1 Publications**

Basic Information

Observed MW

120kDa

Calculated MW

96kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Involved in histone H4-K16 acetylation; histone H4-K5 acetylation; and histone H4-K8 acetylation. Located in nucleoplasm. Part of histone acetyltransferase complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

55683

Swiss Prot

Q9P2N6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human KANSL3 (NP_001108488.1).

Synonyms

NSL3; Rcd1; KIAA1310; KANSL3

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using KANSL3 antibody (A8234) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 15s.