

A8521

Leader in Biomolecular Solutions for Life Science



ANKRD44 Rabbit pAb

Catalog No.: A8521

Basic Information

Observed MW

108kDa

Calculated MW

108kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Putative regulatory subunit of protein phosphatase 6 (PP6 that may be involved in the recognition of phosphoprotein substrates.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

91526

Swiss Prot

Q8N8A2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 310-550 of human ANKRD44 (NP_001182073.1).

Synonyms

PP6-ARS-B; ANKRD44

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using ANKRD44 Rabbit pAb (A8521) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.