

A8580

Leader in Biomolecular Solutions for Life Science



NPHP3 Rabbit pAb

Catalog No.: A8580

Basic Information

Observed MW

151kDa

Calculated MW

151kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

This gene encodes a protein containing a coiled-coil (CC) domain, a tubulin-tyrosine ligase (TTL) domain, and a tetratricopeptide repeat (TPR) domain. The encoded protein interacts with nephrocystin, it is required for normal ciliary development, and it functions in renal tubular development. Mutations in this gene are associated with nephronophthisis type 3, and also with renal-hepatic-pancreatic dysplasia, and Meckel syndrome type 7. Naturally occurring read-through transcripts exist between this gene and the downstream ACAD11 (acyl-CoA dehydrogenase family, member 11) gene.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

27031

Swiss Prot

Q7Z494

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human NPHP3 (NP_694972.3).

Synonyms

MKS7; NPH3; RHPD; RHPD1; SLSN3; CFAP31; NPHP3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using NPHP3 antibody (A8580) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.