

A8623

Leader in Biomolecular Solutions for Life Science



CPLX1 Rabbit pAb

Catalog No.: A8623

Basic Information

Observed MW

Refer to figures

Calculated MW

15kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Mouse

Background

Proteins encoded by the complexin/synaphin gene family are cytosolic proteins that function in synaptic vesicle exocytosis. These proteins bind syntaxin, part of the SNAP receptor. The protein product of this gene binds to the SNAP receptor complex and disrupts it, allowing transmitter release.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

10815

Swiss Prot

O14810

Immunogen

A synthetic Peptide of human CPLX1

Synonyms

CPX1; CPX-I; DEE63; EIEE63; CPLX1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.