

A8657

Leader in Biomolecular Solutions for Life Science



## Septin 9 Rabbit pAb

Catalog No.: A8657

1 Publications

### Basic Information

#### Observed MW

75kDa

#### Calculated MW

65kDa

#### Category

Polyclonal Antibody

#### Applications

WB, IF/ICC, ELISA

#### Cross-Reactivity

Human, Mouse, Rat

### Background

This gene is a member of the septin family involved in cytokinesis and cell cycle control. This gene is a candidate for the ovarian tumor suppressor gene. Mutations in this gene cause hereditary neuralgic amyotrophy, also known as neuritis with brachial predilection. A chromosomal translocation involving this gene on chromosome 17 and the MLL gene on chromosome 11 results in acute myelomonocytic leukemia. Multiple alternatively spliced transcript variants encoding different isoforms have been described.

### Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200

### Immunogen Information

#### Gene ID

10801

#### Swiss Prot

Q9UHD8

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human Septin 9 (NP\_001106964.1).

#### Synonyms

MSF; MSF1; NAPB; SEPT9; SINT1; PNUTL4; SeptD1; AF17q25; Septin 9

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

## Validation Data



Western blot analysis of extracts of various cell lines, using Septin 9 antibody (A8657) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 5s.



Confocal immunofluorescence analysis of U2OS cells using Septin 9 Polyclonal Antibody (A8657) at dilution of 1:100. Blue: DAPI for nuclear staining.