

A8741

Leader in Biomolecular Solutions for Life Science



# RRAS Rabbit pAb

Catalog No.: A8741

## Basic Information

### Observed MW

23kDa

### Calculated MW

23kDa

### Category

Polyclonal Antibody

### Applications

WB,IP

### Cross-Reactivity

Human,Mouse,Rat

## Background

The protein encoded by this gene is a small GTPase involved in diverse processes including angiogenesis, vascular homeostasis and regeneration, cell adhesion, and neuronal axon guidance. Mutations in this gene are found in many invasive cancers.

## Recommended Dilutions

WB 1:500 - 1:2000

IP 1:20 - 1:50

## Immunogen Information

### Gene ID

6237

### Swiss Prot

P10301

### Immunogen

Recombinant protein of human RRAS

### Synonyms

R-Ras; RRAS

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of extracts of human spleen, using RRAS antibody (A8741).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.