

A8897

Leader in Biomolecular Solutions for Life Science



RAB27B Mouse mAb

Catalog No.: A8897 **1 Publications**

Basic Information

Observed MW

27kDa

Calculated MW

25kDa

Category

Monoclonal Antibody

Applications

WB,FC,ELISA

Cross-Reactivity

Human,Mouse

Background

Enables guanyl ribonucleotide binding activity; myosin V binding activity; and protein domain specific binding activity. Involved in multivesicular body sorting pathway and positive regulation of exocytosis. Located in Golgi stack and cytoplasmic vesicle.

Recommended Dilutions

WB 1:500 - 1:2000

FC 1:20 - 1:50

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

5874

Swiss Prot

O00194

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

C25KG; RAB27B

Contact

 www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,pH7.3.