

A8955

Leader in Biomolecular Solutions for Life Science



ACAP2 Rabbit pAb

Catalog No.: A8955

Basic Information

Observed MW

Refer to figures

Calculated MW

88kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human,Mouse

Background

Enables GTPase activator activity. Acts upstream of or within actin filament-based process. Located in ruffle.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

23527

Swiss Prot

Q15057

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

CENTB2; CNT-B2; ACAP2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.