

A8960

Leader in Biomolecular Solutions for Life Science



## FLAD1 Rabbit pAb

Catalog No.: A8960

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

65kDa

#### Category

Polyclonal Antibody

#### Applications

WB

#### Cross-Reactivity

Human

### Background

This gene encodes the enzyme that catalyzes adenylation of flavin mononucleotide (FMN) to form flavin adenine dinucleotide (FAD) coenzyme. Alternatively spliced transcript variants encoding distinct isoforms have been observed.

### Recommended Dilutions

WB 1:500 - 1:1000

### Immunogen Information

#### Gene ID

80308

#### Swiss Prot

Q8NFF5

#### Immunogen

A synthetic Peptide of human FLAD1

#### Synonyms

FAD1; FADS; PP591; LSMFLAD; FLAD1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide, pH7.3.