

A9027

Leader in Biomolecular Solutions for Life Science



LIN7C Rabbit pAb

Catalog No.: A9027

Basic Information

Observed MW

Refer to figures

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables L27 domain binding activity and cytoskeletal protein binding activity. Involved in morphogenesis of an epithelial sheet. Located in cell-cell junction; cytoplasm; and plasma membrane. Part of MPP7-DLG1-LIN7 complex.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

55327

Swiss Prot

Q9NUP9

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

MALS3; VELI3; LIN-7C; MALS-3; LIN-7-C; LIN7C

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.