

A9030

Leader in Biomolecular Solutions for Life Science



AGAP2 Rabbit pAb

Catalog No.: A9030

Basic Information

Observed MW

Refer to figures

Calculated MW

125kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

The protein encoded by this gene belongs to the centaurin gamma-like family. It mediates anti-apoptotic effects of nerve growth factor by activating nuclear phosphoinositide 3-kinase. It is overexpressed in cancer cells, and promotes cancer cell invasion. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

116986

Swiss Prot

Q99490

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

PIKE; GGAP2; CENTG1; AGAP2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.