

A9063

Leader in Biomolecular Solutions for Life Science



# PODN Rabbit pAb

Catalog No.: A9063

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

69kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

The protein encoded by this gene is a member of the small leucine-rich repeat protein family and contains an amino terminal CX3CXCX7C cysteine-rich cluster followed by a leucine-rich repeat domain. Studies suggest that this protein could function to inhibit smooth muscle cell proliferation and migration following arterial injury.

## Recommended Dilutions

WB 1:500 - 1:1000

## Immunogen Information

### Gene ID

127435

### Swiss Prot

Q7Z5L7

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

PCAN; SLRR5A; PODN

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.