

A9063

Leader in Biomolecular Solutions for Life Science



PODN Rabbit pAb

Catalog No.: A9063

Basic Information

Observed MW

Refer to figures

Calculated MW

69kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

The protein encoded by this gene is a member of the small leucine-rich repeat protein family and contains an amino terminal CX3CXCX7C cysteine-rich cluster followed by a leucine-rich repeat domain. Studies suggest that this protein could function to inhibit smooth muscle cell proliferation and migration following arterial injury.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

127435

Swiss Prot

Q7Z5L7

Immunogen

A synthetic Peptide of human PODN

Synonyms

PCAN; SLRR5A; PODN

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.