

A9140

Leader in Biomolecular Solutions for Life Science



KIF2C/MCAK Rabbit mAb

Catalog No.: A9140

Recombinant

Basic Information

Observed MW

85kDa

Calculated MW

81kDa

Category

SMab Recombinant Monoclonal
Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

CloneNo number

ARC1447

Background

This gene encodes a kinesin-like protein that functions as a microtubule-dependent molecular motor. The encoded protein can depolymerize microtubules at the plus end, thereby promoting mitotic chromosome segregation. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

11004

Swiss Prot

Q99661

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIF2C/KIF2C/MCAK (Q99661).

Synonyms

MCAK; CT139; KNSL6; KIF2C/MCAK

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using KIF2C/KIF2C/MCAK Rabbit mAb (A9140) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 3min.