

A9184

Leader in Biomolecular Solutions for Life Science



RIN1 Rabbit pAb

Catalog No.: A9184

Basic Information

Observed MW

Refer to figures

Calculated MW

84kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Predicted to enable small GTPase binding activity. Predicted to be involved in endocytosis; regulation of catalytic activity; and signal transduction. Predicted to act upstream of or within associative learning; memory; and negative regulation of synaptic plasticity. Located in cytoplasm and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

9610

Swiss Prot

Q13671

Immunogen

A synthetic Peptide of human RIN1

Synonyms

RIN1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.