

A9184

Leader in Biomolecular Solutions for Life Science



RIN1 Rabbit pAb

Catalog No.: A9184

Basic Information

Observed MW

Refer to figures

Calculated MW

84kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Predicted to enable small GTPase binding activity. Predicted to be involved in endocytosis; regulation of catalytic activity; and signal transduction. Predicted to act upstream of or within associative learning; memory; and negative regulation of synaptic plasticity. Located in cytoplasm and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

9610

Swiss Prot

Q13671

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

RIN1

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.