

A9457

Leader in Biomolecular Solutions for Life Science



TACO1 Rabbit pAb

Catalog No.: A9457

Basic Information

Observed MW

32kDa

Calculated MW

32kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a mitochondrial protein that function as a translational activator of mitochondrially-encoded cytochrome c oxidase 1. Mutations in this gene are associated with Leigh syndrome.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

51204

Swiss Prot

Q9BSH4

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

CCDC44; MC4DN8; TACO1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

Validation Data



Western blot analysis of lysates from 293T cells, using TACO1 Rabbit pAb (A9457).
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.