

A9535

Leader in Biomolecular Solutions for Life Science



HES5 Rabbit pAb

Catalog No.: A9535

Basic Information

Observed MW

Refer to figures

Calculated MW

18kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

This gene encodes a member of a family of basic helix-loop-helix transcriptional repressors. The protein product of this gene, which is activated downstream of the Notch pathway, regulates cell differentiation in multiple tissues. Disruptions in the normal expression of this gene have been associated with developmental diseases and cancer.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

388585

Swiss Prot

Q5TA89

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

bHLHb38; HES5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.