

A9718

Leader in Biomolecular Solutions for Life Science



# TACC3 Rabbit pAb

Catalog No.: A9718

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

90kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

This gene encodes a member of the transforming acidic coiled-coil protein family. The encoded protein is a motor spindle protein that may play a role in stabilization of the mitotic spindle. This protein may also play a role in growth and differentiation of certain cancer cells.

## Recommended Dilutions

WB 1:500 - 1:1000

## Immunogen Information

### Gene ID

10460

### Swiss Prot

Q9Y6A5

### Immunogen

A synthetic Peptide of human TACC3

### Synonyms

ERIC1; Tacc4; ERIC-1; maskin; TACC3

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.