

A9718

Leader in Biomolecular Solutions for Life Science



TACC3 Rabbit pAb

Catalog No.: A9718

Basic Information

Observed MW

Refer to figures

Calculated MW

90kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes a member of the transforming acidic coiled-coil protein family. The encoded protein is a motor spindle protein that may play a role in stabilization of the mitotic spindle. This protein may also play a role in growth and differentiation of certain cancer cells.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

10460

Swiss Prot

Q9Y6A5

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

ERIC1; Tacc4; ERIC-1; maskin; TACC3

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH 7.3.