

A9984

Leader in Biomolecular Solutions for Life Science



MZT2B Rabbit pAb

Catalog No.: A9984

Basic Information

Observed MW

16kDa

Calculated MW

16kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Located in cytosol; microtubule cytoskeleton; and nucleoplasm. Part of gamma-tubulin large complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

80097

Swiss Prot

Q6NZ67

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human MZT2B (NP_079305.2).

Synonyms

FAM128B; MOZART2B; MZT2B

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using MZT2B antibody (A9984) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.