

A9994

Leader in Biomolecular Solutions for Life Science



# FTSJ3 Rabbit pAb

Catalog No.: A9994

## Basic Information

---

### Observed MW

125kDa

### Calculated MW

97kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human

## Background

---

Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli Ftsj protein and to a mouse protein possibly involved in embryogenesis.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

117246

### Swiss Prot

Q8IY81

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 608-847 of human FTSJ3 (NP\_060117.3).

### Synonyms

SPB1; EPCS3; FTSJ3

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of extracts of various cell lines, using FTSJ3 antibody (A9994) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 30s.