

AC005N

Leader in Biomolecular Solutions for Life Science



Rabbit Control IgG

Catalog No.: AC005N **1 Publications**

Basic Information

Observed MW

Calculated MW

Category

SMab Recombinant Monoclonal
Antibody

Applications

IP

Cross-Reactivity

Background

Recommended Dilutions

IP 1:50 - 1:200

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Contact

 www.abclonal.com

Product Information

Source
Rabbit


Isotype
IgG

Purification
Potein A/G purification

Storage

Store at 4°C Avoid freeze / thaw cycles.
Buffer: PBS,pH7.3.

Validation Data



Immunoprecipitation of over-expressed DDDDK-tagged protein in 293T cells incubated using DDDDK-tag antibody (AE063). A mock served as negative control using rabbit Control IgG (AC005N) and over-expressed 293T cell lysate served as positive control.