

AC034

Leader in Biomolecular Solutions for Life Science



Rat Control IgG

Catalog No.: AC034

Basic Information

Observed MW

Calculated MW

Category

SMab Recombinant Monoclonal
Antibody

Applications

IP

Cross-Reactivity

-

Background

Recommended Dilutions

IP 1:50 - 1:200

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

None

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rat

Isotype

IgG

Purification

Potein A/G purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.