

AE077

Leader in Biomolecular Solutions for Life Science



Rabbit anti GST-Tag mAb

Catalog No.: AE077

6 Publications

Basic Information

Observed MW

27kDa/33.9kDa

Calculated MW

26kDa

CategorySMab Recombinant Monoclonal
Antibody**Applications**

WB,IP,ELISA

Cross-Reactivity

Species independent

CloneNo number

ARC50738

Background

Glutathione S-transferases (GSTs), previously known as ligandins, comprise a family of eukaryotic and prokaryotic phase II metabolic isozymes best known for their ability to catalyze the conjugation of the reduced form of glutathione (GSH) to xenobiotic substrates for the purpose of detoxification. The GST family consists of three superfamilies: the cytosolic, mitochondrial, and microsomal—also known as MAPEG—proteins. Members of the GST superfamily are extremely diverse in amino acid sequence, and a large fraction of the sequences deposited in public databases are of unknown function. The Enzyme Function Initiative (EFI) is using GSTs as a model superfamily to identify new GST functions. A GST-tag is often used to separate and purify proteins that contain the GST-fusion protein. The tag is 220 amino acids (roughly 26 kDa) in size, which, compared to tags such as the Myc-tag or the FLAG-tag, is quite large. It can be fused to either the N-terminus or C-terminus of a protein. However, many commercially available sources of GST-tagged plasmids include a thrombin domain for cleavage of the GST tag during protein purification.

Recommended Dilutions

WB 1:2000 - 1:6000**IP** 3µg antibody for 1µg
extracts of recombinant
protein**ELISA** Recommended starting
concentration is 1
µg/mL. Please optimize
the concentration
based on your specific
assay requirements.

Immunogen Information

Gene ID**Swiss Prot**

P08515

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of schja GST protein (P08515).

Synonyms

GST;GST tag;GST-tag

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

Contact

www.abclonal.com

Validation Data



Western blot analysis of lysates from GST protein, using Rabbit anti GST-Tag mAb (AE077) at 1:5000 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 3s.



Immunoprecipitation analysis of 1 µg extracts of GST Protein using 3 µg Rabbit anti GST-Tag antibody (AE077). Western blot was performed from the immunoprecipitate using Rabbit anti GST-Tag antibody (AE077) at a dilution of 1:1000.