

AE124

Leader in Biomolecular Solutions for Life Science



# Magnetic Beads-conjugated Rabbit anti HA-Tag mAb

Catalog No.: AE124

## Basic Information

### Observed MW

37kDa

### Calculated MW

### Category

SMab Recombinant Monoclonal  
Antibody

### Applications

IP

### Cross-Reactivity

Species independent

### CloneNo number

ARC59578-MB

### Conjugate

Magnetic Beads

## Recommended Dilutions

IP                      20µl-40µl Magnetic  
Beads for 100µg-300µg  
extracts of whole cells

## Background

## Immunogen Information

### Gene ID

Swiss Prot

### Immunogen

A synthetic peptide corresponding to HA tag.

### Synonyms

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG1,Kappa


### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.05% proclin300,pH7.3.

# Validation Data



Immunoprecipitation of LDHA-3X HA in 300 µg extracts from 293T cells transfected with LDHA-3X HA using 20 µl Magnetic Beads-conjugated Rabbit anti HA-Tag mAb (AE124). Western blot analysis was performed using HRP-conjugated Mouse anti HA-Tag mAb (AE025) at 1:5000 dilution.