

AP0206

Leader in Biomolecular Solutions for Life Science



Phospho-Paxillin-Y31 Rabbit pAb

Catalog No.: AP0206

Basic Information

Observed MW

Refer to figures

Calculated MW

65kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
5829

Swiss Prot
P49023

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Paxillin; PXN; paxillin; Phospho-Paxillin-Y31

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.