

AP0350

Leader in Biomolecular Solutions for Life Science



## Phospho-FAM129A-S602 Rabbit pAb

Catalog No.: AP0350

### Basic Information

**Observed MW**

150kDa

**Calculated MW**

103kDa

**Category**

Polyclonal Antibody

**Applications**

WB

**Cross-Reactivity**

Human

### Background

This gene encodes a member of the family with sequence similarity 129 protein family. This gene is highly expressed in several cancer cells and may serve as a prognostic marker for certain cancers. The encoded protein may play a role in regulating p53-mediated apoptosis.

### Recommended Dilutions

WB 1:500 - 1:1000

### Immunogen Information

**Gene ID**

116496

**Swiss Prot**

Q9BZQ8

**Immunogen**

Synthetic peptide. This information is considered to be commercially sensitive.

**Synonyms**

GIG39; NIBAN; C1orf24; FAM129A; Phospho-FAM129A-S602

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.