

AP0376

Leader in Biomolecular Solutions for Life Science



Phospho-JunB-S79 Rabbit pAb

Catalog No.: AP0376

Basic Information

Observed MW

Calculated MW

36kDa

Category

Polyclonal Antibody

Applications

IF/ICC

Cross-Reactivity

Human,Mouse,Rat

Background

Enables sequence-specific double-stranded DNA binding activity. Involved in positive regulation of transcription by RNA polymerase II. Located in nucleoplasm. Part of transcription factor AP-1 complex. Biomarker of Hodgkin's lymphoma and anaplastic large cell lymphoma.

Recommended Dilutions

IF/ICC 1:100 - 1:200

Immunogen Information

Gene ID

3726

Swiss Prot

P17275

Immunogen

A phospho specific peptide corresponding to residues surrounding S79 of human Phospho-JunB-S79

Synonyms

AP-1; Phospho-JunB-S79

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

