

AP0412

Leader in Biomolecular Solutions for Life Science



Phospho-MYL12A-S18 Rabbit pAb

Catalog No.: AP0412 **2 Publications**

Basic Information

Observed MW

20kDa

Calculated MW

20kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human,Rat

Background

This gene encodes a nonsarcomeric myosin regulatory light chain. This protein is activated by phosphorylation and regulates smooth muscle and non-muscle cell contraction. This protein may also be involved in DNA damage repair by sequestering the transcriptional regulator apoptosis-antagonizing transcription factor (AATF)/Che-1 which functions as a repressor of p53-driven apoptosis. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 8.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10627

Swiss Prot

P19105

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

MLCB; MRCL3; MRLC3; MYL2B; MLC-2B; HEL-S-24; Phospho-MYL12A-S18

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

