

AP0602

Leader in Biomolecular Solutions for Life Science



# Phospho-STK11-S428 Rabbit pAb

Catalog No.: AP0602

5 Publications

## Basic Information

### Observed MW

### Calculated MW

49kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse,Rat

## Background

The protein encoded by this gene is a serine/threonine kinase that regulates cell polarity and energy metabolism and functions as a tumor suppressor. Mutations in this gene have been associated with the autosomal dominant Peutz-Jeghers syndrome, as well as with skin, pancreatic, and testicular cancers.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

6794

### Swiss Prot

Q15831

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

PJS; LKB1; hLKB1; Phospho-STK11-S428

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.