

AP0773

Leader in Biomolecular Solutions for Life Science



Phospho-GAB1-Y627 Rabbit pAb

Catalog No.: AP0773

Basic Information

Observed MW

Refer to figures

Calculated MW

77kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
2549

Swiss Prot
Q13480

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

DFNB26; Phospho-GAB1-Y627

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.